

# Sign-in


For existing users, click “Sign In”, enter username and password, and click “Sign In.”



The diagram illustrates the sign-in process. It starts with a 'TeleMeWhere' menu screen featuring a 'Sign In' button (highlighted with a red rectangle) and a 'Create Account' button. An arrow points from the 'Sign In' button to a 'Sign In' form screen. This form screen includes input fields for 'Username' and 'Password', a 'Sign In' button, and a 'Main Menu' link.

# Create User Account

For new users, click “Create Account” to sign up. Enter valid information and then click “Create Account.”






The diagram illustrates the account creation process. It starts with a 'TeleMeWhere' menu screen featuring a 'Sign In' button and a 'Create Account' button (highlighted with a red rectangle). An arrow points from the 'Create Account' button to a 'Create Account' form screen. This form screen includes input fields for 'First Name', 'Last Name', 'Username', 'Password', 'Re-Enter Password', and 'Email', a 'Create Account' button, and a 'Main Menu' link.

- Usage: Use gaze cursor to select buttons. Use pinching gesture to click buttons. Clicking input fields will bring up keyboard. Clicking buttons on keyboard will populate input fields.
- “Main Menu” can be clicked to return to said menu.

# Model manipulation

- Once logged in, Click on any of these action options to proceed to manipulate the model or markers!

Home



# Model Manipulation

- Click “Drag” to be able to drag the model around. This is the default action.
- Click “Rotate” to be able to rotate the model around.
- Click “Homepage” to go back to the doctor/patient homepage.




- Click “Remove Marker” to remove markers by clicking on them.

- Click “Place Marker” and click the model where you wish to place a marker. This will open a menu to customize the marker. By clicking an existing marker, you will be able to view the message attached to that marker.

# Invitation

- Click on “Invite Doctor(Patient)” to grant other users access to your model.
- Click on “Check Invites” to accept or decline other users’ invitation to view their models.

Home



**Patient Homepage**

Invite Doctor

Check Invites

Edit Information

Sign Out

Go Back

**Invite Patient**

Enter the username of the Patient you wish to invite and click "Send Invite"

Username

Send Invite

Go Back

**Invite Patient**

Enter the username of the Patient you wish to invite and click "Send Invite"

Username

- Click on "Invite Doctor(Patient)" to grant other users access to your model.

**Pending Invites**

Invited By...

- Click on "Edit Information" to edit your personal information.

**Patient Homepage**

- Click on "Check Invites" to see what patients or doctors have invited you to view their model.

**Edit Information**

- Click on "To Model" to view and manipulate your own model.